doserec

ARCHIVING AND RETRIEVAL DOSIMETRIC DATA

DoseRecis is a system for collecting dosimetric data from radiological examinations, featuring statistical functions in accordance with the requirements of the Ministerial Decree of 3 November 2023.

DoseRec interfaces via DICOM with the radiological examination archives within healthcare facilities (PACS).

Data can be fed into the system automatically, via pre-set forwards from the PACS, or via the Query/Retrieve module when automatic transmission from the PACS is not possible.

The **DoseRec** application is web-based and accessible via a browser using credentials to protect data.

DATA COLLECTION

In order to store dosimetric data in **DoseRec**, diagnostic departments must send the following to the PACS system:

- the SRDose file, certified as compliant, for CT and mammography examinations
- DICOM images with DAP data saved in dose-specific DICOM tags for X-rays and orthopantomograms

Each examination is recorded in the system, collecting crucial data, including:

- Data on the equipment used, including name and model
- Study ID, description and date
- Radiation dose data (for example, average glandular dose per MG, dose-length product for CT, DAP value for X-rays and orthopantomography)
- Patient data, such as gender, date of birth and ORTHOPANTO (optional in the DICOM file)
- Information on the target anatomical region

DATA PROCESSING

The data is processed to extract meaningful information regarding the radiation doses to which patients were exposed during the examinations. The following aspects are considered:

- **Age Classification:** The patient's age is calculated based on their date of birth and the date of the study. It is then classified into age groups (e.g., 'Paediatric', 'Adult').
- **Classification by Side and Anatomical Region:** Mammograms are analysed based on side (e.g., "Left breast", "Right breast"), whilst CT scans, X-rays and orthopantomograms are classified according to anatomical region (e.g., "Chest", "Abdomen").
- **Determination of Subcategories:** Target anatomical regions are categorised into subcategories using the SNOMED system, to ensure that each piece of data is correctly associated with the relevant anatomical part.

CALCULATION OF STATISTICS AND DISPLAY OF RESULTS

Once the data has been filtered, the following statistics are calculated for each combination of gender, age group, side (for MG) or anatomical region (for CT), and machine used.

The results are displayed in an easily understandable format, with the data broken down by categories. Valid data and excluded data are shown separately, with a brief description of the reason for exclusion in cases where a study does not meet the requirements. The data is also presented with aggregated information on the total dose administered during each study.

The processed data is thus ready to be reported to the relevant ministerial portal.

The screenshot displays the DoseRec web application interface. The browser address bar shows the URL `172.30.0.28:8000/home/`. The page title is "Home".

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Footer: Copyright © 2023 G. SQUARED S.p.A. - Via S. Maria 101 - 20122 Milano - Italia. Tutti i diritti sono riservati.